



Anti-TMEM200B (aa 72-121) polyclonal antibody (CABT-BL3617)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 72-121 (AGYWPHRAGA PGSRAANASS PQMSELRREG RGGGRAHGPV ERLRLGPVI) of Human TTMB
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Membrane; Multi-pass membrane protein
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [399474](#)

Protein Refseq [NP_001003682](#)

UniProt ID [Q69YZ2](#)
