



Anti-TMEM229A (aa 161-190) polyclonal antibody (CABT-BL3619)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide between amino acids 161-190 from the central region of Human TMEM229A, conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane.
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	730130
Protein Refseq	NP_001129474

UniProt ID

[B2RXF0](#)
