



# Anti-TMEM93 (aa 2-51) polyclonal antibody (CABT-BL3627)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                              |                                                                                                                                                                           |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen</b>             | Synthetic peptide within residues: AAVVAKREGP PFISEAAVRG NAAVLDYCRT SVSALSGATA GILGLTGLYG, corresponding to N terminal amino acids 2-51 of Human Transmembrane protein 93 |
| <b>Isotype</b>               | IgG                                                                                                                                                                       |
| <b>Source/Host</b>           | Rabbit                                                                                                                                                                    |
| <b>Species Reactivity</b>    | Human                                                                                                                                                                     |
| <b>Purification</b>          | Immunogen affinity purified                                                                                                                                               |
| <b>Conjugate</b>             | Unconjugated                                                                                                                                                              |
| <b>Applications</b>          | WB                                                                                                                                                                        |
| <b>Cellular Localization</b> | Membrane; Multi-pass membrane protein                                                                                                                                     |
| <b>Format</b>                | Liquid                                                                                                                                                                    |
| <b>Size</b>                  | 100 µl                                                                                                                                                                    |
| <b>Buffer</b>                | 2% Sucrose, PBS                                                                                                                                                           |
| <b>Preservative</b>          | None                                                                                                                                                                      |
| <b>Storage</b>               | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.                                                                            |

## GENE INFORMATION

Entrez Gene ID [83460](#)

---

Protein Refseq [NP\\_001014764](#)

---

UniProt ID [Q9BV81](#)

---