



Anti-TMSB15B (aa 1-45) KLH) polyclonal antibody (CABT-BL3635)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic human Thymosin b15 (aa 1-45) KLH conjugated(MGDRPDLSEVERFDKSKLKKKTITEVKN TLPSKETIEQEKEFVKRS)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	286527
Protein Refseq	NP_919305
UniProt ID	P0CG34