



Anti-TOMM20L polyclonal antibody (CABT-BL3641)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Conjugate	Unconjugated
Applications	WB, ELISA
Buffer	PBS with 0.09% (W/V) sodium azide
Preservative	0.09% Sodium Azide
Storage	Maintain refrigerated at 2-8 °C for up to 6 months. For long term storage store at -20 °C in small aliquots to prevent freeze-thaw cycles.

GENE INFORMATION

Entrez Gene ID	387990
Protein Refseq	NP_997260
UniProt ID	Q6UXN7