



Anti-TPRN polyclonal antibody (CABT-BL3644)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human C9orf75.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	This locus encodes a sensory epithelial protein. It was defined by linkage analysis in three Pakistani families to lie between D9S1818 (centromeric) and D9SH6 (telomeric). Mutations at this locus have been associated with autosomal recessive deafness.
---------------------	---

GENE INFORMATION

Entrez Gene ID [286262](#)

Protein Refseq [NP_001121700](#)

UniProt ID [Q4KMQ1](#)
