



Anti-TRAPPC6B polyclonal antibody (CABT-BL3650)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide (Human)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Golgi apparatus, cis-Golgi network. Endoplasmic reticulum.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	TRAPPC6B is a component of TRAPP complexes, which are tethering complexes involved in vesicle transport (Kummel et al.,;2005 (PubMed 16025134)).
--------------	--

GENE INFORMATION

Entrez Gene ID [122553](#)

Protein Refseq [NP_001073005](#)

UniProt ID [Q86SZ2](#)
