



# Anti-TRCG1 (N-terminal) polyclonal antibody (CABT-BL3656)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A synthetic peptide from the N terminal region of TRCG1 conjugated to an immunogenic carrier protein.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, IHC-P, IHC-Fr
<b>Format</b>	Liquid
<b>Buffer</b>	Whole serum
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles.

## GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">541610</a>
<b>Protein Refseq</b>	<a href="#">NP_001014420</a>

UniProt ID

[F8VPJ5](#)

---