



# Anti-TRDMT1 (aa 39-53&361-376) polyclonal antibody (CABT-BL3658)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptides ANEVYKHNFPTHLL and PEKTTVKQRYRLLGNS, corresponding to amino acids 39-53 and 361-376 of Mouse Dnmt2.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Mouse, Human
<b>Purification</b>	Protein G purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Cellular Localization</b>	Nucleus.
<b>Format</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	Specifically methylates cytosine 38 in the anticodon loop of tRNA(Asp).
---------------------	---

# GENE INFORMATION

Entrez Gene ID	<a href="#">13434</a>
Protein Refseq	<a href="#">NP_034197</a>
UniProt ID	<a href="#">O55055</a>