



Anti-TRIP12 (aa 75-125) polyclonal antibody (CABT-BL3673)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region between residues 75 and 125 of human TRIP12
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µl
Buffer	0.1% BSA, Tris buffered saline
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction	Probable E3 ubiquitin-protein ligase TRIP12 is an enzyme that in humans is encoded by the TRIP12 gene.
--------------	--

GENE INFORMATION

Entrez Gene ID [9320](#)

Protein Refseq [NP_004229](#)

UniProt ID [Q14669](#)
