



Anti-TRP polyclonal antibody (CABT-BL3675)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide from 4th cytoplasmic loop of drosophila trp conjugated to an immunogenic carrier protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Drosophila
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr
Cellular Localization	Cell Membrane; Multi-pass membrane protein
Format	Liquid
Buffer	PBS, Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	43542
Protein Refseq	NP_476768

UniProt ID

[P19334](#)
