



Anti-TRPM6 polyclonal antibody (CABT-BL3679)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: GCERDKNRSSLERHTRL
Isotype	IgG
Source/Host	Guinea pig
Species Reactivity	Rat
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	IHC-Fr
Cellular Localization	Multi-pass membrane protein.
Format	Liquid
Buffer	PBS, Whole serum
Preservative	0.05% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	TRPM6 is a transient receptor potential ion channel associated with hypomagnesemia with secondary hypocalcemia.
---------------------	---

GENE INFORMATION

Entrez Gene ID	225997
Protein Refseq	NP_700466
UniProt ID	Q8CIR4