



Anti-TRPM6 polyclonal antibody (CABT-BL3680)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide (Human) conjugated to immunogenic carrier protein
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-Fr, IHC-P
Cellular Localization	Multi-pass membrane protein.
Format	Liquid
Buffer	PBS
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	TRPM6 is a transient receptor potential ion channel associated with hypomagnesemia with secondary hypocalcemia.
---------------------	---

GENE INFORMATION

Entrez Gene ID [140803](#)

Protein Refseq [NP_001170781](#)

UniProt ID [Q9BX84](#)
