



Anti-TSPAN33 (aa 108-157) polyclonal antibody (CABT-BL3684)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within the internal sequence amino acids 108-157 (QLAAGILGFV FSDKARGKVS EIINNAIVHY RDDLDLQNL DFGQKKFSCC) of human Tetraspanin-33.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Plays an important role in normal erythropoiesis. It has a role in the differentiation of erythroid progenitors.
---------------------	--

GENE INFORMATION

Entrez Gene ID [340348](#)

Protein Refseq [NP_848657](#)

UniProt ID [Q86UF1](#)
