



Anti-TSPAN33 polyclonal antibody (CABT-BL3685)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human TSPAN33.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Plays an important role in normal erythropoiesis. It has a role in the differentiation of erythroid progenitors.
--------------	--

GENE INFORMATION

Entrez Gene ID [340348](#)

Protein Refseq [NP_848657](#)

UniProt ID [Q86UF1](#)
