



Anti-TTLL3 (aa 1-102) polyclonal antibody (CABT-BL3696)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	TEGDRNIWIV KPGAKSRGRG IMCMDHLEEM LKLVNGNPVV MKDGKWWVQK, corresponding to amino acids 1-102 of Human TTLL3
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, ICC/IF
Format	Liquid
Size	50 µl
Buffer	1X PBS, pH 7.2
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [26140](#)

Protein Refseq [NP_001021100](#)

UniProt ID [J3KQB2](#)
