



Anti-TUBGCP6 (internal region) polyclonal antibody (CABT-BL3701)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to internal of human TUBGCP6.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB, IHC, ELISA
Format	Liquid
Buffer	In PBS, pH 7.4 (150mM NaCl, 0.02% sodium azide, 50% glycerol)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The protein encoded by this gene is part of a large multisubunit complex required for microtubule nucleation at the centrosome.
--------------	---

GENE INFORMATION

Entrez Gene ID	85378
----------------	-----------------------

Protein Refseq [NP_065194](#)

UniProt ID [Q96RT7](#)
