



# Anti-TUSC5 (internal region) polyclonal antibody (CABT-BL3703)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A synthetic peptide derived from an internal region of human TUSC5
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA, ICC/IF
<b>Cellular Localization</b>	Cell Membrane; Multipass membrane protein
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	50% Glycerol, PBS, 150mM Sodium chloride, pH 7.4
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at -20°C. Stable for 12 months at -20°C

## GENE INFORMATION

Entrez Gene ID [286753](#)

Protein Refseq [NP\\_758955](#)

---

UniProt ID [Q8IXB3](#)

---