



Anti-TUSC5 (C-terminal) polyclonal antibody (CABT-BL3704)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to C terminal residues of Mouse TUSC5
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Cell Membrane; Multipass membrane protein
Format	Liquid
Buffer	50% Glycerol, PBS
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	237858
Protein Refseq	NP_808377

UniProt ID

[Q8C838](#)
