



Anti-TYW5 (aa 266-315) polyclonal antibody (CABT-BL3716)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 266-315 (NKDPTAASRA AQILDRALKT LAELPEEYRD FYARRMVLHI QDKAYSKNSE) of human C2orf60.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	129450
Protein Refseq	NP_001034782

UniProt ID

[A2RUC4](#)
