



Anti-UBAC2 (aa 108-157) polyclonal antibody (CABT-BL3717)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide to a region within residues: CSIPRVQVAQ ILGPLSITNK TLIYILGLQL FTSGSYIWIV AISGLMSGLC, corresponding to internal sequence amino acids 108-157 of Human UBAC2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Multi-pass membrane protein
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Entrez Gene ID [337867](#)

Protein Refseq [NP_808882](#)

UniProt ID [A0A024RE02](#)
