



# Anti-UBP64E (aa 475-574) polyclonal antibody (CABT-BL3726)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                              |                                                                                                                                                                                                        |
|------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen</b>             | Fusion protein: GRHITVTDLALPDLYKPRVYFYNP SLKKMKITRVYVSQSF<br>NINLV LMSAYEMLNVEQFAPLSRCRLVAYNSSMDTIIQSL ESCTDPALTELRAANYSLDF,<br>corresponding to amino acids 475/574 of Drosophila melanogaster Ubp64E |
| <b>Isotype</b>               | IgG                                                                                                                                                                                                    |
| <b>Source/Host</b>           | Mouse                                                                                                                                                                                                  |
| <b>Species Reactivity</b>    | Fruit fly                                                                                                                                                                                              |
| <b>Purification</b>          | Whole antiserum                                                                                                                                                                                        |
| <b>Conjugate</b>             | Unconjugated                                                                                                                                                                                           |
| <b>Applications</b>          | WB                                                                                                                                                                                                     |
| <b>Cellular Localization</b> | Nuclear                                                                                                                                                                                                |
| <b>Format</b>                | Liquid                                                                                                                                                                                                 |
| <b>Concentration</b>         | Lot specific                                                                                                                                                                                           |
| <b>Buffer</b>                | 50% Glycerol                                                                                                                                                                                           |
| <b>Preservative</b>          | None                                                                                                                                                                                                   |
| <b>Storage</b>               | Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.                                                                                       |

## GENE INFORMATION

Entrez Gene ID [38644](#)

---

Protein Refseq [NP\\_729087](#)

---

UniProt ID [Q24574](#)

---