



Anti-UBQLN4 (internal region) polyclonal antibody (CABT-BL3727)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from internal sequence of Human UBQLN4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P
Cellular Localization	Nucleus. Cytoplasm. Endoplasmic reticulum. May associate with the endoplasmic reticulum.
Format	Liquid
Buffer	50% Glycerol, PBS (without Mg ²⁺ and Ca ²⁺), 150mM Sodium chloride, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	UBQLN4 (ubiquilin 4) is a protein-coding gene. Diseases associated with UBQLN4 include mumps. GO annotations related to this gene include protein binding and polyubiquitin binding. An important paralog of this gene is UBQLN3.
---------------------	---

GENE INFORMATION

Entrez Gene ID [56893](#)

Protein Refseq [NP_064516](#)

UniProt ID [Q9NRR5](#)
