



Anti-UBQLN4 (N-terminal) polyclonal antibody (CABT-BL3728)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to N-terminus 29 amino acids of human UBQLN4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Conjugate	Unconjugated
Applications	WB, IHC, ELISA
Format	Liquid
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for three month. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	UBQLN4 (ubiquilin 4) is a protein-coding gene. Diseases associated with UBQLN4 include mumps. GO annotations related to this gene include protein binding and polyubiquitin binding. An important paralog of this gene is UBQLN3.
--------------	---

GENE INFORMATION

Entrez Gene ID [56893](#)

Protein Refseq [NP_064516](#)

UniProt ID [Q9NRR5](#)
