



Anti-UBTFL1 polyclonal antibody (CABT-BL3729)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to a region of Human UBTFL1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	UBTFL1 is a preimplantation-specific gene and is involved in early development, implantation, and embryonic stem (ES);cell derivation (summary by Yamada et al., 2010 (PubMed 19915186)).
---------------------	---

GENE INFORMATION

Entrez Gene ID [642623](#)

Protein Refseq [NP_001137447](#)

UniProt ID [P0CB47](#)
