



# Anti-UGT2B11 (aa 428-458) polyclonal antibody (CABT-BL3732)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 428-458 of Human UGT2B11.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Cellular Localization</b>	Microsome membrane. Endoplasmic reticulum membrane.
<b>Format</b>	Liquid
<b>Buffer</b>	99% PBS
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at 4°C (up to 6 months). For long term storage store at -20°C

## GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">10720</a>
<b>Protein Refseq</b>	<a href="#">NP_001064</a>

UniProt ID

[O75310](#)

---