



# Anti-UMODL1 (aa 1115-1164) polyclonal antibody (CABT-BL3734)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                              |                                                                                                                                                                       |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen</b>             | Synthetic peptide corresponding to a region within internal sequence amino acids 1115-1164 (EMQLFIGDSP IPQNYSVSAS DDVRIEVGLY RQKSNLKVV L TECWATPSSN) of Human UMODL1. |
| <b>Isotype</b>               | IgG                                                                                                                                                                   |
| <b>Source/Host</b>           | Rabbit                                                                                                                                                                |
| <b>Species Reactivity</b>    | Human                                                                                                                                                                 |
| <b>Purification</b>          | Immunogen affinity purified                                                                                                                                           |
| <b>Conjugate</b>             | Unconjugated                                                                                                                                                          |
| <b>Applications</b>          | WB                                                                                                                                                                    |
| <b>Cellular Localization</b> | Cytoplasm and Cell membrane.                                                                                                                                          |
| <b>Format</b>                | Liquid                                                                                                                                                                |
| <b>Buffer</b>                | 2% Sucrose, PBS                                                                                                                                                       |
| <b>Preservative</b>          | None                                                                                                                                                                  |
| <b>Storage</b>               | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.                                                                        |

## GENE INFORMATION

Entrez Gene ID [89766](#)

Protein Refseq [NP\\_001004416](#)

---

UniProt ID [Q5DID0](#)

---