



Anti-UNK polyclonal antibody (CABT-BL3739)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human UNK.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	UNK unkempt is the human homolog of the Drosophila melanogaster RING finger protein unkempt. This transcription factor contains multiple zinc fingers and coiled-coil domains. In Drosophila, it is essential during post-embryonic development. Its function in mammalian development is unclear.
---------------------	--

GENE INFORMATION

Entrez Gene ID [85451](#)

Protein Refseq [NP_001073888](#)

UniProt ID [A0A024R8N4](#)
