



# Anti-URM1 (N-terminal) polyclonal antibody (CABT-BL3743)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A synthetic peptide (conjugated with KLH) corresponding to N-terminus of human URM1.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Protein A purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, IHC
<b>Format</b>	Liquid
<b>Buffer</b>	In PBS (0.09% sodium azide)
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	Acts as a sulfur carrier required for 2-thiolation of mcm(5)S(2)U at tRNA wobble positions of tRNA(Lys), tRNA(Glu) and tRNA(Gln). Serves as sulfur donor in tRNA 2-thiolation reaction by thiocarboxylated (-COSH) at its C-terminus by MOCS3. The sulfur is t
---------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

# GENE INFORMATION

Entrez Gene ID [81605](#)

---

Protein Refseq [NP\\_001129419](#)

---

UniProt ID [Q9BTM9](#)

---