



Anti-UTP18 (aa 325-375) polyclonal antibody (CABT-BL3752)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region between residue 325 and 375 of human UTP18.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IP
Cellular Localization	Nucleus, nucleolus.
Format	Liquid
Buffer	8mM PBS, 60mM Citrate, 150mM Tris, pH 7-8
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction	U3 small nucleolar RNA-associated protein 18 homolog is a protein that in humans is encoded by the UTP18 gene.
--------------	--

GENE INFORMATION

Entrez Gene ID [51096](#)

Protein Refseq [NP_057085](#)

UniProt ID [Q9Y5J1](#)
