



Anti-VMA21 (aa 1-50) polyclonal antibody (CABT-BL3765)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A region within synthetic peptide: MERPDKAALN ALQPPEFRNE SSLASTLKTL LFFTALMITV PIGLYFTTKS, corresponding to N terminal amino acids 1-50 of Human VMA21
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P
Cellular Localization	Endoplasmic reticulum
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	This gene encodes a vacuolar ATPase assembly integral membrane protein. Mutations of this gene result in an X-linked vacuolar myopathy with excessive autophagy.
---------------------	--

GENE INFORMATION

Entrez Gene ID [203547](#)

Protein Refseq [NP_001017980](#)

UniProt ID [Q3ZAQ7](#)
