



# Anti-VMN1R43 polyclonal antibody (CABT-BL3769)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A synthetic peptide from mouse vomeronasal type-1 receptor A5 (VN1A5, V1ra5, V1ra2, V1RA2) conjugated to an immunogenic carrier protein was used as the antigen.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Purification</b>	Whole serum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC, WB
<b>Format</b>	Lyophilised
<b>Preservative</b>	None
<b>Storage</b>	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.
<b>Ship</b>	This item will be shipped to you at ambient temperature in a lyophilised form.

## GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">113847</a>
<b>Protein Refseq</b>	<a href="#">NP_444450.2</a>

