



Anti-VMN1R44 polyclonal antibody (CABT-BL3770)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide as a part of mouse vomeronasal type-1 receptor B4 (VN1B4, V1ra10, V1rb11) and vomeronasal type-1 receptor B2 (VN1B2, V1rb2, VRi2) conjugated to immunogenic carrier protein was used as the immunogen.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Whole serum
Conjugate	Unconjugated
Applications	IHC, WB
Format	Lyophilised
Preservative	None
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.
Ship	This item will be shipped to you at ambient temperature in a lyophilised form.

GENE INFORMATION

Entrez Gene ID [113854](#)

Protein Refseq [NP_444457.2](#)

UniProt ID Q9EQ47
