



Anti-VMN1R46 polyclonal antibody (CABT-BL3772)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|------------------------------|--|
| Immunogen | Synthetic peptide as a part of mouse vomeronasal type-1 receptor B8 or A13 (VN1B8, V1rb8, V1ra13) conjugated to immunogenic carrier protein. |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Mouse |
| Purification | Whole antiserum |
| Conjugate | Unconjugated |
| Applications | WB, IHC-P, IHC-Fr |
| Cellular Localization | Cell Membrane; Multi-pass membrane protein |
| Format | Liquid |
| Buffer | Whole serum |
| Preservative | None |
| Storage | Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles. |

GENE INFORMATION

Entrez Gene ID [113856](#)

Protein Refseq [NP_444459](#)

UniProt ID [Q9EQ45](#)
