



# Anti-VMN1R50 polyclonal antibody (CABT-BL3775)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A synthetic peptide from mouse vomeronasal type-1 receptor B1 (VN1B1, V1rb1, V1ra5, VN2, V1RA5) conjugated to an immunogenic carrier protein was used as the immunogen.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Purification</b>	Whole serum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC, WB
<b>Format</b>	Lyophilised
<b>Preservative</b>	None
<b>Storage</b>	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.
<b>Ship</b>	This item will be shipped to you at ambient temperature in a lyophilised form.

## GENE INFORMATION

Entrez Gene ID [113852](#)

Protein Refseq [NP\\_444455.1](#)

