



# Anti-VMN1R53 polyclonal antibody (CABT-BL3777)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                              |  |
|------------------------------|--|
| <b>Immunogen</b>             | Synthetic peptide as a part of mouse vomeronasal type-1 receptor B3 (VN1B3, V1rb3, VN5, pheromone receptor VN5) conjugated to immunogenic carrier protein          |
| <b>Isotype</b>               | IgG  |
| <b>Source/Host</b>           | Rabbit   |
| <b>Species Reactivity</b>    | Mouse  |
| <b>Purification</b>          | Whole antiserum  |
| <b>Conjugate</b>             | Unconjugated   |
| <b>Applications</b>          | ELISA, IHC-P, IHC-Fr   |
| <b>Cellular Localization</b> | Cell Membrane  |
| <b>Format</b>                | Liquid   |
| <b>Buffer</b>                | Whole serum  |
| <b>Preservative</b>          | None   |
| <b>Storage</b>               | Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles. |

## GENE INFORMATION

Entrez Gene ID [113853](#)

Protein Refseq [NP\\_444456](#)

---

UniProt ID [A2RS76](#)

---