



Anti-VOM1R105 polyclonal antibody (CABT-BL3780)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide as a part of rat vomeronasal type-1 receptor A1, A13, A14 (VN1A1, Vnr1, V1A13, V1ra13, V1A14, V1ra14) conjugated to immunogenic carrier protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr
Cellular Localization	Cell Membrane; Multi-pass membrane protein
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction putative pheromone receptor; may play a role in olfactory sensory perception.

GENE INFORMATION

Entrez Gene ID [286893](#)

Protein Refseq [NP_775136](#)

UniProt ID [Q62850](#)
