



Anti-VOM1R94 polyclonal antibody (CABT-BL3782)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide as a part of rat vomeronasal type-1 receptor A1, A13,B25(VN1A1, Vnr1, V1A13, V1ra13, V1A14, V1ra14) conjugated to immunogenic carrier protein was used as the immunogen.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Purification	Whole serum
Conjugate	Unconjugated
Applications	IHC, WB
Format	Lyophilised
Preservative	None
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.
Ship	This item will be shipped to you at ambient temperature in a lyophilised form.

GENE INFORMATION

Entrez Gene ID [494265](#)

Protein Refseq [NP_001008915.1](#)

UniProt ID Q5J3M4
