



Anti-VSTM4 (aa 149-177) polyclonal antibody (CABT-BL3785)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, between amino acids 149-177 from the Central region of Human C10orf72 isoform 1 and 2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell membrane.
Format	Liquid
Size	80 µl
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID [196740](#)

Protein Refseq [NP_659421](#)

UniProt ID [Q8IW00](#)
