



Anti-W02A11.3 (aa 28-127) polyclonal antibody (CABT-BL3788)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: CMLESLIRRPFCMNPQLQTHPYDRQPTTAPRLQMMPRRDQ TQNPVIAELMSTVTRVRDRKRQRDEEAMYPRRSVGLQTD EDVSEVLDGGTRISDMKRVR, corresponding to amino acids 28/127 of C. elegans 1N331
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Caenorhabditis elegans
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	173144
Protein Refseq	NP_493231

UniProt ID

[Q9XUM8](#)
