



Anti-WDR70 polyclonal antibody (CABT-BL3798)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region between residue 604 and 654 of Human WDR70
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Buffer	0.1% BSA, Tris buffered saline
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

GENE INFORMATION

Entrez Gene ID	55100
Protein Refseq	NP_060504
UniProt ID	Q9NW82