



Anti-WDR8 (aa 146-195) polyclonal antibody (CABT-BL3801)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 146-195 (RYLALAERRD CRDYVSIFVC SDWQLLRHFD TDTQDLTGIE WAPNGCVLAA) of Mouse WDR8.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasm.
Format	Liquid
Buffer	97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction WD repeat-containing protein 8 is a protein that in humans is encoded by the WDR8 gene.

GENE INFORMATION

Entrez Gene ID [59002](#)

Protein Refseq [NP_067474](#)

UniProt ID [Q9JM98](#)
