



Anti-wrm-1 (aa 6-105) polyclonal antibody (CABT-BL3806)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: AETFSQPCTPLNFPMTPTSRVSTPVRPSSTMSARQYSG SPFKAQPQNMEPSNSRVQELREAAVGKRSYTNAWMQGTYA PPQMAGQQRFPSRPPSVIGS, corresponding to amino acids 6/105 of Caenorhabditis elegans Wrm1
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Caenorhabditis elegans
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [175802](#)

Protein Refseq [NP_498236](#)

UniProt ID [Q10953](#)
