



# Anti-WUN2 (aa 1-90) polyclonal antibody (CABT-BL3807)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                              |                                                                                                                                                                                                          |
|------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen</b>             | Fusion protein: MSTLRPVSVCTDTPLQRFESQSSSGEEPSSPTASSIIAAA SGPNNNNNNN VKLDLQLPTFVDSTRSQGSL YAVASNKKPP LRGPRQIFGR, corresponding to amino acids 1/90 of Fruit fly ( <i>Drosophila melanogaster</i> ) Wunen2 |
| <b>Isotype</b>               | IgG                                                                                                                                                                                                      |
| <b>Source/Host</b>           | Mouse                                                                                                                                                                                                    |
| <b>Species Reactivity</b>    | Fruit fly                                                                                                                                                                                                |
| <b>Purification</b>          | Whole antiserum                                                                                                                                                                                          |
| <b>Conjugate</b>             | Unconjugated                                                                                                                                                                                             |
| <b>Applications</b>          | WB                                                                                                                                                                                                       |
| <b>Cellular Localization</b> | Plasma membrane                                                                                                                                                                                          |
| <b>Format</b>                | Liquid                                                                                                                                                                                                   |
| <b>Buffer</b>                | 50% Glycerol                                                                                                                                                                                             |
| <b>Preservative</b>          | None                                                                                                                                                                                                     |
| <b>Storage</b>               | Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.                                                                                         |

## GENE INFORMATION

Entrez Gene ID [53558](#)

Protein Refseq [NP\\_524991](#)

---

UniProt ID [Q7KGH1](#)

---