



Anti-XIRP1 (aa 1337-1367) polyclonal antibody (CABT-BL3810)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 1337-1367 of Human XIRP1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell junction. Colocalizes with actin stress fibers.
Format	Liquid
Buffer	99% PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	165904
Protein Refseq	NP_919269

UniProt ID

[Q702N8](#)
