



Anti-XRRA1 (aa 1-50) polyclonal antibody (CABT-BL3820)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to N terminal amino acids 1 - 50 (MAFSGIYKLD DGKPYLNNCF PARNLLRVPE EGQGHWLVVQ KGNLKKKPKG) of Human XRRA1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic and Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	143570
Protein Refseq	NP_892014

UniProt ID

[Q6P2D8](#)
