



Anti-YIPF4 (aa 35-84) polyclonal antibody (CABT-BL3824)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 35-84 (PDVKLNLLGGD FIKESTATT F LRQRGYGWLL EVEDDDPEDN KPLLEELDID) of Human YIPF4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane.
Format	Liquid
Buffer	98% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Protein YIPF4 is a protein that in humans is encoded by the YIPF4 gene.
---------------------	---

GENE INFORMATION

Entrez Gene ID [84272](#)

Protein Refseq [NP_115688](#)

UniProt ID [Q9BSR8](#)
