



Anti-YIPF6 (aa 179-228) polyclonal antibody (CABT-BL3825)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: NFMVRLFVVI VMFAWSIVAS TALLADSQPP NRRALAVYPV FLFYFVISWM, corresponding to C terminal amino acids 179-228 of Human YIPF6
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane; Multi-pass membrane protein
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Protein YIPF6 is a protein that in humans is encoded by the YIPF6 gene.
---------------------	---

GENE INFORMATION

Entrez Gene ID [286451](#)

Protein Refseq [NP_001182143](#)

UniProt ID [Q96EC8](#)
