



Anti-ZACN (aa 36-85) polyclonal antibody (CABT-BL3827)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 36 - 85 (SLFNVNLSKK VQESIQIPNN GSAPLLVDVR VFVSNVFNVD ILRYTMSSML) of Human LGICZ1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Integral to plasma membrane; postsynaptic membrane.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	LGICZ1 is a zinc-activated ligand-gated ion channel that defines a new subgroup of the cysteine-loop superfamily of ligand-gated ion channels (Davies et al., 2003 (PubMed
---------------------	--

12381728)).

GENE INFORMATION

Entrez Gene ID	353174
----------------	------------------------

Protein Refseq	NP_851321
----------------	---------------------------

UniProt ID	Q401N2
------------	------------------------
