



Anti-ZBTB10 (aa 50-100) polyclonal antibody (CABT-BL3828)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region within amino acids 50-100 of Human ZBT10.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IP
Cellular Localization	Nucleus.
Format	Liquid
Size	100 µg
Buffer	99% Tris citrate/phosphateNote: pH 7-8
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction	May be involved in transcriptional regulation.
---------------------	--

GENE INFORMATION

Entrez Gene ID [65986](#)

Protein Refseq [NP_001099009](#)

UniProt ID [Q96DT7](#)
